



PRICE LIST

EFFECTIVE JANUARY 1st, 2008

(ITEMS & PRICES SUBJECT TO CHANGE WITHOUT NOTICE OR OBLIGATION)



102-122 WEST MAIN STREET • LISBON, IA 52253-0899

319-455-2110 • 800-553-7297 • FAX 319-455-2166

E-mail: sales@lloydtable.com • Web: www.lloydtable.com

	<u>Page</u>		<u>Page</u>
Galaxy 906HS, 900HS, Hylo, Elevation, Stationary,	2	Cosmic III Intersegmental Traction Table	9
Options for Galaxy Tables	3	Astro Tables & Options	10 & 11
Ultimate, McManis Elevation, Hylo & Stationary	4	Side Posture Tables	11
Options for McManis Tables	5	ACCENTS / Activator Tables	12
402 Flexion Elevation / Stationary & Options	6 & 7	Adjusting Benches, Pediatric, Portables	13 & 14
Upper Cervical Elevation Table	7	Chairs, Knee Chest, SOT, Exam Tables	15
400 Series & 402 Elevation/Stationary & Options	8 & 9	General Information, Guarantee	16

GALAXY 906 H. S.

BASE PRICE: **10695⁰⁰**



Shown With Optional Drops

TABLE SPECIFICATIONS:

Height Range	18" to 38"
Length - Extended	82 Inches
Length - Retracted	70 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	450 Lbs.
Table Wt (approx)	400 Lbs.

STANDARD FEATURES:

Deluxe/Adjustable Headpiece
 One Piece Adjustable Chest With Adjustable
 Tension Spring Breakaway
 Chest/Pelvic Separation
 Extension/Retraction of Length (Power)
 Adjustable Pelvic (Vertical)(Power)
 Adjustable Ankle Rest (Vertical)
 Speed Adjustment of Hylo & Elevation
 Electronic Height Selector
 Able to Tilt From Any Height In As Little
 As 9 Seconds

GALAXY 900 H. S.

BASE PRICE: **10450⁰⁰**



Shown With Optional Drops

TABLE SPECIFICATIONS:

Height Range	19" to 36"
Length - Extended	84 Inches
Length - Retracted	72 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	460 Lbs.
Table Wt (approx)	410 Lbs.

STANDARD FEATURES:

Deluxe/Adjustable Headpiece
 One Piece Adjustable Chest With Adjustable
 Tension Spring Breakaway
 Chest/Pelvic Separation
 Extension/Retraction of Length (Power)
 Adjustable Pelvic (Vertical)
 Adjustable Ankle Rest (Vertical)
 Speed Adjustment of Hylo Feature
 Tilt Safety Switch
 Long Lasting Naugahyde Upholstery (see pg.16)

GALAXY HYLO

BASE PRICE: **8025⁰⁰**



TABLE SPECIFICATIONS:

Height Range	18" to 28"
Length - Extended	84 Inches
Length - Retracted	72 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	320 Lbs.
Table Wt (approx)	270 Lbs.

STANDARD FEATURES:

Deluxe/Adjustable Headpiece
 One Piece Adjustable Chest With Adjustable
 Tension Spring Breakaway
 Chest/Pelvic Separation
 Extension/Retraction of Length (Power)
 Adjustable Pelvic (Vertical)
 Adjustable Ankle Rest (Vertical & Horizontal)
 Speed Adjustment of Hylo Feature
 Long Lasting Naugahyde Upholstery (see pg.16)

GALAXY ELEVATION / STATIONARY

BASE PRICE: **6925⁰⁰/3025⁰⁰**



Shown With Optional Drops

TABLE SPECIFICATIONS:

Height Range-Elevation	19" to 36"
Height Range-Elevation 06	17" to 37"
Height Range-Stationary	18" to 28"
Length - Extended	83 Inches
Length - Retracted	62 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	400/200 Lbs.

STANDARD FEATURES:

Deluxe/Adjustable Headpiece
 One Piece Adjustable Chest With Adjustable
 Tension Spring Breakaway
 Chest/Pelvic Separation (06 - Powered)
 Extension/Retraction of Length (Manual)
 Adjustable Pelvic (Vertical) (06 - Powered)
 Adjustable Ankle Rest (Vertical & Horizontal)
 Hydraulic Height Variation (Elevation Only)

GALAXY ELEVATION 06 (Not Pictured)

7925⁰⁰

GUARANTEE

LLOYD equipment is guaranteed against defects in parts, materials (except upholstery) and workmanship as follows:

Galaxy Tables, Astro Tables, Activator Tables, 402 Tables, Cosmic III:

ONE YEAR (Full Parts)

SECOND & THIRD YEAR (Limited Parts)

LIFETIME

THREE YEARS

Adjusting Bench Frames

All Other Products:

Guarantee void if warranty card not returned to Lloyd Table Company

OPTIONAL EQUIPMENT FOR GALAXY TABLES

HEAD SECTION:

Straight Drop Headpiece	\$ 450.00
Adjustable Forward Motion & Straight Drop Head (Manual / Automatic)	\$ 770.00 / 1000.00
CBP Type Head (Translation, Forward & Reverse Drop (Manual / Auto??).....	\$ 1370.00 / 1600.00
7" Elevation, Tilt, Lateral Bending, Adj. Width Cushions with Built-in Paper)	
Total Dimension Flexion™	2000.00
Includes Forward & Straight Drop, Elevation, Rotation, Lateral Bending, Manual Flexion, Manual Extension & Retraction, Cervical Flexion Assist Handle & Harness	

CHEST/LUMBAR SECTION:(All-New T-Style / Total Response Drop System Now Standard)

Chest and Lumbar Drop(One Piece Cushion) (Manual)	\$ 1240.00
Chest and Lumbar Drop (One Piece Cushion)(Manual & Automatic)(<u>Auto Drops Standard W/O Selector Valve</u>)	1620.00
Pre 2005 Drop Style Chest & Lumbar(2-Piece Cushions), ADD to Drop Charge	245.00
Spring Breakaway Chest With Adjustable Tension (Add if Chest & Lumbar Drop is Ordered)	375.00
Chest Rotation (Includes Spring Breakaway-without Adjustable Tension):	
When Standard Chest is Ordered	1215.00
When Chest/Lumbar Drop is Ordered (Add to Drop Charge)	945.00
CBP Type Chest (Add to Drop Charge)	655.00
Standard Chest (2 Cushions - No Drops)	95.00
Split Chest with Filler	110.00
Breast Cushion (Center Section Firm Foam, Outer Sections Softer Foam)	40.00

PELVIC SECTION & OTHER OPTIONS:(All-New T-Style / Total Response Drop System Now Standard)

Pelvic Drop (Manual)	\$ 1050.00
Pelvic Drop (Manual & Automatic):	
Galaxy 900 HS, Galaxy Hylo, Galaxy Elevation	1975.00
Multi-Directional Pelvic Drop (Add to Manual or Manual & Automatic Drop Charge)	740.00
Cephalic Drops of 52°, 39°, 26° & 13°; Straight Down Drop, Caudal Drops of 13°, 26°, 39° & 52°	
Drop Ankle (Manual / Automatic).....	\$ 715.00 / 1265.00
Electronic Height Selector - Galaxy Elevation / 06	\$ 975.00
Semi-Automatic Foot Pedal - Galaxy Elevation / 06.....	395.00
Fold-in Footplate - Galaxy 900 HS, Galaxy Hylo	215.00

REPLACEMENT STAPLE-ON AND/OR SNAP-ON COVERS

Complete Set Staple-On and/or Snap-On Covers:		
With Standard Chest/One Piece Drop Chest.....	\$ 215.00 / 235.00	Headpiece (Pair).....
With Chest/Lumbar Drop(2-Piece Cushion).....	230.00 / 250.00	Standard Chest
Cloth Upholstery - Table or Snap-On Covers - Add.....	75.00	Chest Drop.....
		Lumbar Drop
		Pelvic.....
		Armrest (Pair).....
		Ankle Rest.....

WIDE (21") CUSHIONS ARE STANDARD - CALL FOR PRICING ON NON-STANDARD WIDTHS

WHEN ORDERING: SPECIFY HEIGHT & INDICATE WITH OR WITHOUT GONADAL NOTCH (No Additional Charge)

AUTO DROP / SELECTOR NOTES

3-Way Selector Valve (Selections for Auto Cocking include:Cervical Only / All Sections / Chest, Lumbar, Pelvic)	\$135.00
8 Position Selector Valve (Pre 2005 Automatic Selector Valve Cocking System)	435.00
If only Auto Head and/or Auto Chest / Lumbar & NO Auto Drop Pelvic is Ordered, Add.....	285.00

GALAXY MODEL TABLES OPTIONAL TOPS

ACTIVATOR TOP	<u>BASE PRICE</u>	203 TYPE TOP	<u>BASE PRICE</u>
47" FACE OPENING	900 HS	18" FACE OPENING	900 HS.....
24"WIDE	HYLO	PAPER ATTACHMENT	HY.....
	ELEVATION	21" WIDE	ELEVATION.....
	\$8750.00		\$ 8750.00
	6000.00		6000.00
	6000.00		6000.00

Call for Additional Cushion Configurations and Dimensions

Items & Prices Subject to Change Without Notice Or Obligation

GALAXY ULTIMATE

Elevation/Hylo/Flexion/Distractation



Shown With Optional Drops

TABLE SPECIFICATIONS:

Height Range	22" to 38"
Length - Extended	91 Inches
Length - Retracted	
(Including Footplate)	81 Inches
(Excluding Footplate)	65 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	500 Lbs.
Table Wt (approx)	450 Lbs.

STANDARD FEATURES:

Electronic Height Selector & Auto Flexion
Deluxe/Adjustable Headpiece
Lumbar Extension Cushion

BASE PRICE:

12995⁰⁰

STANDARD FEATURES CONT:

One Piece Adjustable Chest With Adjustable Tension Spring Breakaway
Chest/Pelvic Separation (Power)
Extension/Retraction of Length (Power)
Adjustable Ankle Rest (Vertical & Horizontal)
Pelvic Section:
 Flexion Distractation (Man.,Semi-Auto, Auto)
 Lateral Bending
 Axial Extension/Retraction (Power)
Ankle Straps
Speed Adjustment of Hylo Feature
Tilt Safety Switch
Long Lasting Naugahyde Upholstery (see pg.16)

GALAXY McMANIS ELEVATION



Shown With Optional Drops

TABLE SPECIFICATIONS:

Height Range	22" to 38"
	24" to 40"
Length - Extended	90 Inches
Length - Retracted	64 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	430 Lbs.
Table Wt (approx)	380 Lbs.

BASE PRICE:

9725⁰⁰

STANDARD FEATURES:

Deluxe/Adjustable Headpiece
Lumbar Extension Cushion
One Piece Adjustable Chest With Adjustable Tension Spring Breakaway
Chest/Pelvic Separation
Extension/Retraction of Length
Pelvic Section:
 Flexion Distractation (Manual & Semi-Auto)
 Lateral Bending
 Axial Extension/Retraction (Power)
Adjustable Ankle Rest (Vertical & Horizontal)
Ankle Straps
Hydraulic/Semi-Automatic Height Variation
Long Lasting Naugahyde Upholstery (see pg.16)

GALAXY McMANIS HYLO



TABLE SPECIFICATIONS:

Heights Available	22" to 28"
Length - Extended	91 Inches
Length - Retracted	
(Including Footplate)	81 Inches
(Excluding Footplate)	65 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	400 Lbs.
Table Wt (approx)	350 Lbs.

BASE PRICE:

9995⁰⁰

STANDARD FEATURES CONT:

STANDARD FEATURES:

Deluxe/Adjustable Headpiece
One Piece Adjustable Chest With Adjustable Tension Spring Breakaway

Lumbar Extension Cushion
Chest/Pelvic Separation
Extension/Retraction of Length (Power)
Pelvic Section:
 Flexion Distractation (Manual & Semi-Auto)
 Lateral Bending
 Axial Extension/Retraction (Power)
Adjustable Ankle Rest (Vertical & Horizontal)
Ankle Straps
Automatic Drop Away Footplate
Speed Adjustment of Hylo Feature
Tilt Safety Switch
Long Lasting Naugahyde Upholstery (see pg.16)

GALAXY McMANIS STATIONARY



TABLE SPECIFICATIONS:

Heights Available	18" to 30"
Length - Extended	90 Inches
Length - Retracted	64 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	220 Lbs.
Table Wt (approx)	170 Lbs.

BASE PRICE:

5125⁰⁰

STANDARD FEATURES CONT:

STANDARD FEATURES:

Deluxe/Adjustable Headpiece
Lumbar Extension Cushion

One Piece Adjustable Chest With Adjustable Tension Spring Breakaway
Chest/Pelvic Separation
Extension/Retraction of Length (Manual)
Pelvic Section:
 Flexion Distractation (Manual)
 Lateral Bending
 Axial Extension/Retraction (Man. Foot Pump)
Adjustable Ankle Rest (Vertical & Horizontal)
Ankle Straps
Long Lasting Naugahyde Upholstery (see pg.16)

OPTIONAL EQUIPMENT FOR GALAXY McMANIS TABLES

HEAD SECTION:

Straight Drop Headpiece	\$ 450.00
Adjustable Forward Motion & Straight Drop Head (Manual / Automatic)	\$ 770.00 / 1000.00
CBP Type Head (Translation, Forward & Reverse Drop (Manual / Auto??).....	\$ 1370.00 / 1600.00
7" Elevation, Tilt, Lateral Bending, Adj. Width Cushions with Built-in Paper)	
Total Dimension Flexion™	2000.00
Includes Forward & Straight Drop, Elevation, Rotation, Lateral Bending, Manual Flexion & Manual Extension & Retraction, Cervical Flexion Assist Handle & Harness	

CHEST/LUMBAR SECTION:(All-New T-Style / Total Response Drop System Now Standard)

Chest and Lumbar Drop(One Piece Cushion) (Manual)	\$ 1240.00
Chest and Lumbar Drop(One Piece Cushion)(Manual & Automatic)(Auto Drops Standard W/O Selector Valve)	1620.00
Pre 2005 Drop Style Chest & Lumbar (2-Piece Cushions), ADD to Drop Charge	245.00
Spring Breakaway Chest With Adjustable Tension (Add if Chest/Lumbar Drop is Ordered)	375.00
Chest Rotation (Includes Spring Breakaway-without Adjustable Tension):	
When Standard Chest is Ordered	1215.00
When Chest/Lumbar Drop is Ordered (Add to Drop Charge)	945.00
CBP Type Chest (Add to Drop Charge)	655.00
Standard Chest (2 Cushions - No Drops)	95.00
Split Chest with Filler	110.00
Breast Cushion (Center Section Firm Foam, Outer Sections Softer Foam)	40.00

PELVIC SECTION & OTHER OPTIONS: (All-New T-Style / Total Response Drop System Now Standard)

Pelvic Drop (Manual)	\$ 1125.00
Pelvic Drop (Manual & Automatic):	
Ultimate, McManis Elevation, McManis Hylo	2050.00
Pelvic Rotation	355.00
Drop Ankle (Manual Only)	\$ 715.00
Electronic Height Selector - McManis Elevation	\$ 630.00
Automatic Flexion - McManis Elevation & McManis Hylo	900.00
Automatic Flexion & Height Select - McManis Elevation	1400.00
(Includes Auto Pause - Delay of 0 to 30 Seconds in Upper & Lower Limits of Auto Pelvic Flexion)	
Telescoping T-Bar	180.00
Angled T-Bar	250.00
Grab Bar	275.00
Replacement Ankle Straps (Pair)	80.00

REPLACEMENT STAPLE-ON AND/OR SNAP-ON COVERS	
Headpiece	\$43.00 / 46.00
Complete Set Staple-On and/or Snap-On Covers:	
Standard Chest.....	50.00 / 55.00
With Standard Chest/One Piece Drop Chest	\$230.00 / 250.00
With Chest/Lumbar Drop (2-Piece Cushion).....	250.00 / 265.00
Cloth Upholstery - Table or Snap-On Covers - Add	75.00
Lumbar Extension.....	31.00 / 33.00
Chest Drop.....	36.00 / 39.00
Lumbar Drop.....	40.00 / 42.00
Pelvic.....	53.00 / 57.00
Ankle Rest.....	31.00 / 33.00
Armrest (Pair).....	41.00 / 44.00

WIDE (21") CUSHIONS STANDARD - CALL FOR PRICING ON NON-STANDARD WIDTHS

WHEN ORDERING: SPECIFY HEIGHT AND INDICATE WITH OR WITHOUT GONADAL NOTCH (No Additional Charge)

AUTO DROP / SELECTOR NOTES

3-Way Selector Valve (Selections for Auto Cocking include:Cervical Only / All Sections / Chest, Lumbar, Pelvic)	\$135.00
8 Position Selector Valve (Pre 2005 Automatic Selector Valve Cocking System)	435.00
If only Auto Head and/or Auto Chest / Lumbar & NO Auto Drop Pelvic is Ordered, Add.....	285.00

All Items & Prices Subject To Change Without Notice Or Obligation

402 FLEXION ELEVATION

BASE PRICE: **5150⁰⁰**



TABLE SPECIFICATIONS:

Height Range	19.5" to 39"
Length - Extended	73 Inches
Length - Retracted	62 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	300 Lbs.
Table Wt (approx)	240 Lbs.

STANDARD FEATURES:

- Tilt Head w/ Add-On Paper Attachment
- Adjustable Chest
- Pelvic Section:
 - Flexion Distraction (Manual)
 - Lateral Bending
 - Axial Extension/Retraction (Powered)
- Adjustable Ankle Rest (Horizontal)
- Ankle Straps
- Long Lasting Naugahyde Upholstery (see pg.16)

203 Type Top (One Piece Head & Chest) Base Price \$4900.00

402 FLEXION STATIONARY

BASE PRICE: **3900⁰⁰**



TABLE SPECIFICATIONS:

Heights Available	18" to 30"
Length - Extended	78 Inches
Length - Retracted	64 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	220 Lbs.
Table Wt (approx)	170 Lbs.

STANDARD FEATURES:

- Tilt Head w/ Add-On Paper Attachment
- Adjustable Chest
- Pelvic Section:
 - Flexion Distraction (Manual)
 - Lateral Bending
 - Axial Extension/Retraction (Man. Foot Pump)
- Adjustable Ankle Rest (Horizontal)
- Ankle Straps
- Long Lasting Naugahyde Upholstery (see pg.16)

OPTIONAL EQUIPMENT FOR 402 FLEXION TABLES

HEAD SECTION:

402 Flexion Elevation and Stationary

Adjustable Width, Tilting Head with Add-on Paper Attachment	190.00
Adjustable Width, Tilting Head with Built-in Paper Attachment*BIPA	300.00
Adjustable Width, Tilting, Forward Motion with *BIPA	630.00
Adjustable Width, Tilting, Elevating with *BIPA	725.00
Adjustable Width, Tilting, Elevating, Straight Drop Headpiece with *BIPA	840.00
Adjustable Width, Tilting, Elevating, Forward Motion & Straight Drop Head with *BIPA (Manual)	965.00
Manual & Automatic (402 Flexion Elevation Only-See Auto Drop Notes Page 7)	1200.00
CBP Type Head (Translation, Forward & Reverse Drop (Manual / Auto??).....\$	1725.00 / 1975.00
7" Elevation, Tilt, Lateral Bending, Adj. Width Cushions with Built-in Paper)	

Total Dimension Flexion™	2200.00
Includes Forward & Straight Drop, Elevation, Rotation, Lateral Bending, Manual Flexion,	
Manual Extension & Retraction, Cervical Flexion Assist & Harness	
Retrofitable to Lloyd Tables With 1989 & Newer Headpieces - Call With Serial Number To Verify Compatibility & Pricing	

CHEST/LUMBAR SECTION:(All-New T-Style / Total Response Drop System Now Standard)

Split Chest with Filler	100.00
Breast Cushion (Center Section Firm Foam, Outer Sections Softer Foam)	40.00
Chest & Lumbar Drop (One Piece Cushion) (Manual)	1240.00
Chest & Lumbar Drop (One Piece Cushion)(Man & Automatic) (Auto Drops Standard W/O Selector Valve)	1620.00
Pre 2005 Drop Style Chest & Lumbar (2-Piece Cushions), ADD to Drop Charge	245.00
Spring Breakaway Chest/Adjustable Tension (With Standard Chest or Drop Chest & Lumbar	375.00
Chest Rotation (Includes Spring Breakaway-Adjustable Tension):	
With Standard Chest	1215.00
With Chest/Lumbar Drops (Add to Drop Charge)	945.00
CBP Type Chest (Add to Drop Charge)	655.00

402 FLEXION SERIES TABLE OPTIONS CONTINUE ON PAGE 7

OPTIONAL EQUIPMENT FOR 402 FLEXION SERIES TABLES

PELVIC SECTION & OTHER OPTIONS: (All-New T-Style / Total Response Drop System Standard)

Pelvic Drop:	
Manual	\$1125.00
Manual & Automatic (402 Flexion Elevation Only)	2050.00
Pelvic Rotation (Not Available When Automatic Axial Distraction is Ordered)	355.00
Vertical Adjustable Ankle Rest	190.00
Drop Ankle(Manual Only).....	715.00
Automatic Flexion (402 Flexion Elevation Only)	1700.00
Automatic Flexion & EHA (402 Flexion Elevation Only)	2250.00
EHA (Electronic Height Adjustment) (402 Flexion Elevation Only)	550.00
Automatic Axial Distraction (With Auto Flexion)/(Without Auto.Flexion).....	550.00/900.00
Telescoping T-Bar	180.00
Angled T-Bar.....	250.00
Grab Bar	275.00
Replacement Ankle Straps	80.00

REPLACEMENT STAPLE-ON AND/OR SNAP-ON COVERS

Complete Set Staple-On and/or Snap-On Covers:	
With Standard Chest/One Piece Drop Chest.....	\$215.00 / 235.00
With Chest/Lumbar Drop(2-Piece Cushions).....	\$230.00 / 250.00
Cloth Upholstery - Table or Snap-On Covers - Add.....	75.00
Headpiece	\$ 43.00 / 46.00 Pelvic.....53.00 / 57.00
Standard Chest.....	50.00 / 55.00 Ankle Rest.....31.00 / 33.00
Chest Drop.....	36.00 / 39.00 Armrest (Pair).....41.00 / 44.00
Lumbar Drop.....	40.00 / 42.00

WIDE (21") CUSHIONS STANDARD - CALL FOR PRICING ON NON-STANDARD WIDTHS

WHEN ORDERING: SPECIFY HEIGHT AND INDICATE WITH OR WITHOUT GONADAL NOTCH (No Additional Charge)

AUTO DROP / SELECTOR NOTES

3-Way Selector Valve (Selections for Auto Cocking include:Cervical Only / All Sections / Chest, Lumbar, Pelvic) \$	135.00
8 Position Selector Valve (Pre 2005 Automatic Selector Valve Cocking System)	435.00
If only Auto Head and/or Auto Chest / Lumbar & NO Auto Drop Pelvic is Ordered, Add.....	285.00

FOAM: A variety of foam densities are available for doctor/patient preference.

STANDARD FOAM DENSITIES: (As Used In The Past) Soft (35 Lb.); Medium (45 Lb.); Medium Firm (50 Lb.); Firm (70 Lb.); Extra Firm (90 Lb. Urethane)

SPECIAL ORDER FOAM DENSITIES: XXX Firm (Bonded); Maximum Firm (Closed Cell Foam) - *Call for Pricing - No Exchanges or Returns on Special Order Foam.*

UPPER CERVICAL SPECIFIC ELEVATION BASE PRICE: **2895⁰⁰**


TABLE SPECIFICATIONS:

Height	9" to 27"
Length	63 Inches
Width	20 Inches
Ship Wt (approx)	180 Lbs.
Table Wt (approx)	120 Lbs.

STANDARD FEATURES:

- Specifically Designed Headpiece with
- Easily Adjustable Tilt and Elevation
- Gray Baked Urethane Powder Coated Frame
- Long Lasting Naugahyde Upholstery (see pg.16)

OPTIONS: Electronic Height Adjustment.....	\$ 550.00
Hideaway Adjustable Ankle Rest (Horizontal)	150.00
Upper Cervical Specific Stationary w/Platform.....	1750.00
Upper Cervical Specific Portable.....	895.00



400 SERIES BENCHES & TABLES

402 ELEVATION-No Drops (Call for Specifications)

BASE PRICE: **3950⁰⁰**

402 ELEVATION-Full Drops

BASE PRICE: **4950⁰⁰**



TABLE SPECIFICATIONS:

Height Range	19" to 38"
Foam	3" Firm
Length - Extended	75 Inches
Length - Retracted	62 Inches
Table Width	22 Inches
Shipping Wt (approx)	280 Lbs.
Table Wt (approx)	225 Lbs.

STANDARD FEATURES: (With Drops)

Adjustable Width Tilting Headpiece
 Adjustable Chest & Pelvic to Accommodate Pregnant & Large Abdominal Patients
 Newly Designed *Total Drop System*TM
 Forward Motion Drop Head w/ Built-In Paper Chest Drop, Hand Cocking
 Pelvic/Lumbar Drop, Hand Cocking
 Adjustable Ankle Rest (Horizontal)
 Long Lasting Naugahyde Upholstery (see pg.16)

402 STATIONARY-No Drops (Call for Specifications)

BASE PRICE: **1850⁰⁰**

402 STATIONARY-Full Drops

BASE PRICE: **2850⁰⁰**



TABLE SPECIFICATIONS:

Adjustable Hts. Available	18" to 26"
Foam	3" Firm
Length	70 Inches
Table Width	22 Inches
Shipping Wt (approx)	165 Lbs.
Table Wt (approx)	115 Lbs.

New!
OPTIONAL KICK DROPS
 See Pricing on Page 9

STANDARD FEATURES:

Adjustable Width Tilting Headpiece
 Adjustable Chest & Pelvic to Accommodate Pregnant & Large Abdominal Patients
 Newly Designed *Total Drop System*TM
 Forward Motion Drop Head w/ Built-In Paper Chest Drop, Hand Cocking
 Pelvic/Lumbar Drop, Hand Cocking
 Adjustable Ankle Rest (Horizontal)
 Adjustable Height Legs
 Long Lasting Naugahyde Upholstery (see pg.16)

401 DB (Drop Bench)

BASE PRICE: **1850⁰⁰**



*Shown with optional kick drops

TABLE SPECIFICATIONS:

Adjustable Hts. Available	18" to 26"
Foam	3" Firm
Length-Extended	72 Inches
Length-Retracted	61 Inches
Table Width	22 Inches
Shipping Wt (approx)	150 Lbs.
Table Wt (approx)	130 Lbs.

STANDARD FEATURES:

Adjustable Width Tilting Headpiece
 Newly Designed *Total Drop System*TM
 Forward Motion Drop Head w/ Built-In Paper Chest Drop, Hand Cocking
 Pelvic/Lumbar Drop, Hand Cocking
 Adjustable Ankle Rest (Horizontal)
 Adjustable Height Legs
 Long Lasting Naugahyde Upholstery (see pg.16)

400 AH

BASE PRICE: **575⁰⁰**



TABLE SPECIFICATIONS:

Adjustable Hts. Available	18" to 26"
Foam	3" Firm
Length - Extended	72 Inches
Length - Retracted	61 Inches
Table Width	22 Inches
Shipping Wt (approx)	150 Lbs.
Table Wt (approx)	130 Lbs.

STANDARD FEATURES:

Tilting Headpiece - Add-on Paper Attachment
 Arm Rests
 Adjustable Height Legs
 Long Lasting Naugahyde Upholstery (see pg.16)

OPTIONAL EQUIPMENT FOR 400 SERIES BENCHES & TABLES

HEAD SECTION:

Adjustable Width, Tilting with Add-on Paper Attachment (400 AH, 402 Stationary & Elevation w/o drops)	\$ 190.00
Adjustable Width, Tilting with Built-in Paper Attachment *BIPA(400 AH, 402 Stat & Elev w/o drops)	300.00
Adjustable Width, Tilting, Forward Motion Drop with *BIPA (400 AH, 402 Stationary & Elevation w/o drops) ...	630.00
Elevating, Adjustable Width, Tilting with *BIPA (400 AH, 402 Stationary & Elevation w/o drops)	725.00
Elevating, Adjustable Width, Tilting, Straight Drop with *BIPA (400 AH, 402 Stationary & Elevation w/o drops) .	840.00
Elevating, Adjustable Width, Tilting, Straight & Forward Motion Drop with *BIPA (400 AH Only)	965.00
Elevating, Straight Drop with *BIPA (401 DB, 402 Stationary & Elevation with drops only)	450.00
CBP Type Head (Translation, Forward & Reverse Drop (Manual / Auto??).....	\$ 1725.00 / 1975.00
7" Elevation, Tilt, Lateral Bending, Adj. Width Cushions with Built-in Paper)	

Total Dimension Flexion™ with Add-on Paper Attachment (400 AH, 401 DB, 402 Stationary & Elevation Only) 2200.00
 Includes Forward & Straight Drop, Elevation, Rotation, Lateral Bending, Manual Flexion & Manual Extension/Retraction, Cervical Flexion Assist Harness

CHEST/LUMBAR/PELVIC SECTIONS:

Drop Pelvic, Hand Cocking (400 AH Only)	\$ 450.00
Drop Chest, Drop Lumbar/Pelvic, Hand Cocking (400AH Only)	770.00
Drop Chest or Drop Lumbar/Pelvic, Kick Drops (400 AH Only) (Add to Cost of Hand Cocking Drops)	ea. 195.00
Kick Drops on 401 DB & 402 Stationary	330.00
Spring Breakaway Chest (402 Stationary & Elevation Only)	300.00
CBP Type Chest (Add to Drop Charge) (402 Stationary & Elevation Only).....	655.00
Auto Cocking (402 Elevation Only)(Includes Auto Chest & Pelvic Only)	950.00

OTHER OPTIONS:

Drop Ankle (Manual / Automatic) (Auto Drop Ankle available 402 Elevation Only).....	\$ 615.00 / 1265.00
EHA (Electronic Height Adjustment) (402 Elevation Only)	550.00
Vertical Adjustable Ankle Rest	190.00
Complete Set Staple-On And/Or Snap-On Covers.....	\$ 215.00 / 235.00

INTERSEGMENTAL TRACTION TABLE

COSMIC III

PRICE:

2250.⁰⁰



Optional 25" Wide Table..... 125.00

TABLE SPECIFICATIONS

Height	25"
Length	72.5"
Width	22"
Table Wt (approx.)	118 lbs.
Shipping Wt (approx.)	175 lbs.

STANDARD FEATURES:

- Clockwise & Counterclockwise Roller Action
- Raise/Lower Rollers up to 4" by 1/10"
- Roller Travel up to 24"
- 0 - 30 Minute Timer
- Automatic Roller Recess
- Vibration
- Digital Electronics
- Hand Control Pad

GUARANTEE

LLOYD equipment is guaranteed against defects in parts, materials (except upholstery) and workmanship as follows:

Galaxy Tables, Astro Tables, Activator Tables, 402 Tables, Cosmic III:	ONE YEAR (Full Parts)
Adjusting Bench Frames	SECOND & THIRD YEAR (Limited Parts)
All Other Products:	LIFETIME
	THREE YEARS

Guarantee void if warranty card not returned to Lloyd Table Company

ASTRO HYLO

BASE PRICE: **5350⁰⁰**



TABLE SPECIFICATIONS:

Heights Available	18" to 28"
Length - Extended	87 Inches
Length - Retracted	76 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	280 Lbs.
Table Wt (approx)	230 Lbs.

STANDARD FEATURES:

- Tilt Headpiece w/Add-on Paper Attachment
- One Piece Chest (Adjustable Tilt)
- Chest/Pelvic Separation
- Ext./Ret. of Length (Manual)
- Adjustable Pelvic (Vertical)
- Adjustable Ankle Rest (Vertical)
- Tilt Safety Switch
- Single Foot Pedal
- Gray Baked Urethane Powder Coated Frame
- Long Lasting Naugahyde Upholstery (see pg.16)

ASTRO ELEVATION

BASE PRICE: **2750⁰⁰**



TABLE SPECIFICATIONS:

Height Range	18" to 38"
Length - Extended	73 Inches
Length - Retracted	60 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	210 Lbs.

STANDARD FEATURES:

- Tilt Headpiece w/Add-on Paper Attachment
- One Piece Body (Chest/Pelvic)
- Adjustable Ankle Rest (Horizontal)
- Gray Baked Urethane Powder Coated Frame
- Long Lasting Naugahyde Upholstery (see pg.16)

202 - Type Top	Base Price	\$2425.00
203 - Type Top	Base Price	\$2525.00

ASTRO STATIONARY KD (KICK DROP)

BASE PRICE: **1695⁰⁰**



Shown With Optional Drops

TABLE SPECIFICATIONS:

Heights Available	18" to 28"
Length - Extended	75 Inches
Length - Retracted	62 Inches
Cushion Width	21 Inches
Table Width	22 Inches
Shipping Wt (approx)	190 Lbs.
Table Wt (approx)	140 Lbs.

STANDARD FEATURES:

- Tilt Headpiece w/Add-on Paper Attachment
- One Piece Body (Chest/Pelvic)
- Adjustable Ankle Rest (Horizontal)
- Long Lasting Naugahyde Upholstery (see pg.16)

OPTIONAL EQUIPMENT FOR ASTRO TABLES

HEAD SECTION:

Astro Hylo, Elevation, and Stationary

Adjustable Width, Tilting Head with Add-on Paper Attachment	190.00
Adjustable Width, Tilting Head with Built-in Paper Attachmen*BIPA	300.00
Adjustable Width, Tilting, Forward Motion with *BIPA	630.00
Adjustable Width, Tilting, Elevating with *BIPA	725.00
Adjustable Width, Tilting, Elevating, Straight Drop Headpiece with *BIPA	840.00
Adjustable Width, Tilting, Elevating, Forward Motion & Straight Drop Head w/ *BIPA (Manual)	965.00
Manual & Automatic Drop Head (Only Available on Asto Hylo when Auto Pelvic is Ordered)	1200.00
CBP Type Head (Translation, Forward & Reverse Drop (Manual / Auto??).....\$	1725.00 / 1975.00
7" Elevation, Tilt, Lateral Bending, Adj. Width Cushions with Built-in Paper)	
Total Dimension Flexion™ with Add-on Paper Attachment	2200.00

(See page 2 for Complete Description)

CHEST/LUMBAR SECTION:

Convert One Piece Body to Chest & Pelvic Cushions	\$ 110.00
Two Cushion Standard Chest - No Drops	105.00
Split Chest with Filler	110.00
Breast Cushion (Center Section Firm Foam, Outer Sections Softer Foam)	40.00
Chest Drop (Manual) (Elevation Only)	410.00
Lumbar Drop (Manual) (Elevation Only)	410.00
Chest & Lumbar Drop New T-Style(One Piece Cushion) (Manual) Astro Stationary Only.....	875.00

OPTIONAL EQUIPMENT FOR ASTRO TABLES

Chest & Lumbar Drop (Manual T-Style - Astro Hylo Only).....	\$1240.00
Chest & Lumbar Drop (Manual & Automatic) (T-Syle Drop W/O Selector Valve) Astro Hylo Only)	1620.00
Pre 2005 Drop Style Chest & Lumbar (Astro Hylo Only)(2-Piece Cushions), ADD to Drop Charge	245.00
Spring Breakaway Chest/Adjustable Tension (With Standard Chest or Drop Chest & Lumbar - Astro Hylo	375.00
Chest Rotation (Includes Spring Breakaway):	
With Standard Chest	1215.00
With Chest/Lumbar Drops (Add to Drop Charge)	945.00
CBP Type Chest (Add to Drop Charge) (Astro Hylo Only).....	655.00

PELVIC SECTION:

Adjustable Pelvic (Astro Elevation & (Astro Stationary w/o a Drop Pelvic) Tables Only) \$ 190.00
Pelvic Drop: Astro Stationary (New T-Style Drop which includes Adjustable Pelvic (Vertical).....	715.00
Astro Elevation (Manual)	525.00
Astro Hylo (T-Style Drop)(Manual / Automatic)	1050.00 / 1975.00

OTHER OPTIONS:

Vertical Adjustable Ankle Rest	\$ 190.00
Drop Ankle (Manual / Automatic) (Manual Only on Astro Stationary).....	\$ 715.00 / 1265.00
EHA (Electronic Height Adjustment) (Astro Elevation & Nicholas Elevation Only)	550.00
Fold-in Footplate (Astro Hylo Only)	215.00
Complete Set Staple-On And/Or Snap-On Covers:	
With Standard Chest.....	\$215.00 / 235.00
With Chest/Lumbar Drop.....	230.00 / 250.00
Cloth Upholstery - Table or Snap-On Covers - Add	75.00

AUTO DROP / SELECTOR NOTES

3-Way Selector Valve (Selections for Auto Cocking include:Cervical Only / All Sections / Chest, Lumbar, Pelvic) \$	135.00
8 Position Selector Valve (Pre 2005 Automatic Selector Valve Cocking System)	435.00
If only Auto Head and/or Auto Chest / Lumbar & NO Auto Drop Pelvic is Ordered, Add.....	285.00

NICHOLAS SIDE POSTURE TABLE

BASE PRICE:

2195⁰⁰




TABLE SPECIFICATIONS:

Height	15 Inches
Length	70 Inches
Width	18 Inches
Ship Wt (approx)	170 Lbs.
Table Wt (approx)	110 Lbs.

STANDARD FEATURES:

Adjustable Headpiece
3 Drop Depths
Sliding Headpiece
Adjustable Shoulder Cushion
Long Lasting Naugahyde Upholstery (see pg.16)

ELEVATING NICHOLAS SIDE POSTURE TABLE

BASE PRICE:

4595⁰⁰



TABLE SPECIFICATIONS:

Height Range	13.5 to 33 Inches
Length	72 Inches
Width	18 Inches
Ship Wt (approx)	210 Lbs.
Table Wt (approx)	260 Lbs.

STANDARD FEATURES:

Hydraulic Elevation
Adjustable Headpiece
3 Drop Depths
Sliding Headpiece
Adjustable Shoulder Cushion
Long Lasting Naugahyde Upholstery (see pg.16)

KALE-TYPE KNEE-CHEST TABLE (Call For Specifications)

BASE PRICE:

1250⁰⁰

GUARANTEE

LLOYD equipment is guaranteed against defects in parts, materials (except upholstery) and workmanship as follows:
Galaxy Tables, Astro Tables, Activator Tables, 402 Tables, Cosmic III:

Adjusting Bench Frames
All Other Products:

ONE YEAR (Full Parts)
SECOND & THIRD YEAR (Limited Parts)
LIFETIME
THREE YEARS

Guarantee void if warranty card not returned to Lloyd Table Company



TABLE SPECIFICATIONS:

Height	24" to 35" Inches
Length	72 Inches
Width	30 Inches
Face Opening	28 Inches
Ship Wt (approx)	37 Lbs.
Table Wt (approx)	31 Lbs.

STANDARD FEATURES:

- Adjustable Height legs
- Seperate Ankle Cushion
- Long Lasting Naugahyde Upholstery (see pg.16)
- UPS Shipping \$ 40.00

INTEGRATOR PORTABLE (Torque Release Technique) BASE PRICE: 565⁰⁰

ACTIVATOR® ADJUSTING BENCH 26" WIDE: 600⁰⁰



TABLE SPECIFICATIONS:

Height	30 Inches
Length	72 Inches
Face Opening	48 Inches
Ship Wt (approx)	150 Lbs.
Table Wt (approx)	130 Lbs.

Other Heights Available Upon Request

STANDARD FEATURES:

- 2" Raised Wedge at Foot End
- 2 ½" Foam
- Long Lasting Naugahyde Upholstery (see pg.16)

Optional 28" WIDE675⁰⁰

ACTIVATOR® ELEVATION BASE PRICE: 3300⁰⁰



TABLE SPECIFICATIONS:

Height	18" to 38"
Length	72 Inches
Width	26 Inches
Face Opening	48 Inches
Ship Wt (approx)	330 Lbs.
Table Wt (approx)	260 Lbs.

STANDARD FEATURES:

- 2" Raised Wedge at Foot End or Separate Ankle Rest
- Gray Baked Urethane Powder Coated Frame

IF ADJUSTABLE ANKLE NOT ORDERED:
SPECIFY RAISED WEDGE AT FOOT END
OR SEPARATE ANKLE REST

OPTIONS: EHA (Electronic Height Adjustment).....\$550.00
 Manual or Automatic Drop Pelvic 1050.00 / 1975.00
 Adjustable Vertical Ankle Rest 280.00

28" Wide Cushions...\$ 125.00

ACTIVATOR® HYLO (PULSAR) BASE PRICE: 5150⁰⁰



TABLE SPECIFICATIONS:

Height	34 Inches
Length	72 Inches
Width	26 Inches
Face Opening	48 Inches
Ship Wt (approx)	240 Lbs.
Table Wt (approx)	195 Lbs.

STANDARD FEATURES:

- Quiet, Smooth Hydraulic Tilting
- Up & Down Tilt Speed Are Individually & Infinitely Adjustable
- Casters for Easy Moving
- Convenient Paper Tray
- Automatic Drop Away Foot Plate
- Gray Baked Urethane Powder Coated Frame
- Long Lasting Naugahyde Upholstery (see pg.16)

OPTIONS: Hand Bars.....\$ 100.00
 28" Wide Cushions..... 125.00
 ATF (Auto Tilt Feature)..... 600.00
 Adjustable Vertical Ankle Rest..... 280.00

IF ADJUSTABLE ANKLE NOT ORDERED:
SPECIFY RAISED WEDGE AT FOOT END
OR SEPARATE ANKLE REST

ACTIVATOR® 906HS (Hilo & Elevation) BASE PRICE: 8995⁰⁰



TABLE SPECIFICATIONS:

Height Range	18" to 38"
Length	84 Inches
Cushion Width	26 Inches
Table Width	25 Inches
Shipping Wt (approx)	425 Lbs.
Table Wt (approx)	375 Lbs.

STANDARD FEATURES:

- Adjustable Vertical Ankle Rest
- Auto Height Select
- Drop-away Foot Plate

OPTIONAL FEATURES:

- Hand Bars..... \$100.00
- 28" Wide Cushions.....125.00
- Manual or Automatic Drop Pelvic... 1050 / 1975

PEDIATRIC PRODUCTS



202 ABSC



203ABNAHC
* w/ Optional
Dinosaur Logo

ADJUSTING BENCHES

202 ABSC - 2" Foam, Face Opening, 48" Long x 15" Wide	Base Price \$ 360.00
203 ABNC - 2" Foam, Face Opening, Narrowed Head & Shoulder Section, 48" Long x 15" Wide	Base Price 470.00
203 ABNAHC - 3" Foam, Face Opening, Narrowed Head & Shoulder Section	
Built-in Paper Attachment, Adjustable Headpiece, 48" Long x 15" Wide	Base Price 635.00
<i>Options for Children's 200 Series Benches - See "Adjusting Bench Options" Above</i>	

401 DBC (Pediatric Drop Bench) BASE PRICE: **1200.⁰⁰**



Height Range	18" to 27"
Length	48"
Width	20"

STANDARD FEATURES:
 Toggle Drop Headpiece
 Chest/Lumbar/Pelvic Drops
 Dinosaur Logo

C 106 Portable - 21" High	260.00
C 106A Portable - 21" to 26"	300.00
Dinosaur Logo For C106 & C106A	45.00
Pediatric Portable Drop Headpiece - 6" x 6" x 2.5"	175.00
Child's Pelvic Adjusting Blocks	
Non-slip or corrugated rubber bottoms	45.00
Pediatric Cushion Top.....	110.00



C106 Child's Portable Table
*Shown with optional Dinosaur Logo

PORTABLE TABLES



C 104

C 105

C 103

C 103 (20½" or 23½" High) (16" Wide)	385.00	UPS	\$20.00
C 104 Tilt Head, Adj. Anklerest, 21" High	450.00	UPS	27.00
C 104A Adjustable Height, 21" to 26"	510.00	UPS	27.00
* Add C 108 Head to C 104 (A).....	395.00	UPS	30.00
C 105 20" Wide, 21" High	395.00	UPS	22.00
C 105A Adjustable Height, 21" to 26"	455.00	UPS	22.00
New! C 107 Tilt Head, Drop Pelvic, 21" High	765.00	UPS	35.00
DNFT Portable	545.00	UPS	27.00
Add-On Paper Attachment for Portables	22.00	UPS	8.00
Zipper Cover for "C" Portable Tables	75.00	UPS	8.00
Integrator Portable	565.00	UPS	40.00

C 108 DROP PORTABLE BASE PRICE: **1420⁰⁰**



TABLE SPECIFICATIONS:

	Standard	Folded
Height	21"	10.5"
Length	67.5"	34"
Width	20"	20"
Ship Wt (approx)	46 Lbs.	
Table Wt (approx)	42 Lbs.	

STANDARD FEATURES:
 Tilting and Elevating Headpiece
 Newly Designed *Full Drop System*™
 Forward Motion & Straight Drop Head
 Chest Drop
 Pelvic/Lumbar Drop

OPTIONS:

Adjustable Height (21 to 26 Inches)	\$ 80.00
Adjustable Ankle Rest (Horizontal)	110.00
With Tilt Head Only (37 Pounds)	1050.00

Items & Prices Subject to Change Without Notice Or Obligation

ADJUSTING BENCHES

BENCH SPECIFICATIONS:	STANDARD FEATURES-ALL BENCHES:
Heights - Adj. by 1" Increments 18" to 28"	Adjustable Height Leg Inserts Are Chrome Plated
Limited Height Adjustment	Gray Baked Urethane Powder Coated Legs
Length 71 Inches	Nauga Plus Upholstery
Width 22 Inches	
Ship Weight 100 Lbs.	

- 201 ABP** Base Price \$ 510.00
 Standard Features: 2" Foam, Face Opening, 2 Cloth Pads, Rectangular Pillow
- 202 AB** Base Price \$ 445.00
 Standard Features: 2" Foam, **NO** Face Opening
- 202 ABS** Base Price \$ 465.00
 Standard Features: 2" Foam, Face Opening
- 203 ABN** Base Price \$ 505.00
 Standard Features: 2" Foam, Face Opening, Narrowed Head & Shoulder Section
- 203 ABN DELUXE**..... Base Price \$705.00
 Standard Features: 3" Foam, Face Opening, Narrowed Head & Shoulder Section,
 Recessed Paper Attachment, Arm Rests, Adjustable Ankle Rest
- 203 ABN AH** Base Price \$ 705.00
 Standard Features: 3" Foam, Face Opening, Narrowed Head & Shoulder Section,
 Recessed Paper Attachment, Adjustable (Tilting) Headpiece
- 203 ABN DELUXE AH** Base Price \$ 800.00
 Standard Features: 3" Foam, Face Opening, Narrowed Head & Shoulder Section,
 Recessed Paper Attachment, Adjustable (Tilting) Headpiece, Arm Rests, Adjustable Ankle Rest

ADJUSTABLE HEIGHT LEGS

Heights Between 20" and 26" With Paper Attachment - May Adjust 2" to 3" Lower or Higher Than Ordered Height
 Heights Between 20" and 26" Without Paper Attachment - May Adjust 4" to 5" Lower or Higher Than Ordered Height
 Heights 18" or Lower - May Not Have Adjustable Heights Available - Call For Information
 *Leg Inserts Available in Lengths of 8" to 18" by 2" Increments

SPECIFY DESIRED HEIGHT WHEN ORDERING

ADJUSTING BENCH OPTIONS

Paper Attachment (Built-in)	\$ 80.00	Heights over 28"	Call for Pricing
Paper Attachment (Add-on)	46.00	Foam Thickness:	
Armrest (Pair)	52.00	2" to 3" or 3" to 4"	\$ 25.00
Adjustable Ankle Rest	100.00	2" to 4"	50.00
Partial Snap-On Cover		Cloth Upholstery	53.00
Body Only 42"	65.00	Adjustable Width Headpiece	
Rectangular Pillow (15" or 18")	38.00	on AH Bench Only	215.00
2 Cloth Pads	35.00		

REPLACEMENTS

Replacement Covers	\$ 135.00
Adjustable Headrest Cover	48.00
Footrest Cover	21.00
Armrest Covers (Pair)	25.00
Armrest Cushions (Pair)	46.00
Adjustable Ht. Legs With Inserts (Pair) ..	220.00
*Adjustable Ht. Leg Inserts (4)	62.00

GUARANTEE

LLOYD equipment is guaranteed against defects in parts, materials (except upholstery) and workmanship as follows:
 Galaxy Tables, Astro Tables, Activator Tables, 402 Tables, Cosmic III:

Adjusting Bench Frames
 All Other Products:

ONE YEAR (Full Parts)
SECOND & THIRD YEAR (Limited Parts)
LIFETIME
THREE YEARS

Guarantee void if warranty card not returned to Lloyd Table Company

Items & Prices Subject To Change Without Notice Or Obligation

CHAIRS • STOOLS • KNEE CHEST • SPINALYZER

201 ABP (See Benches for Details)	\$ 510.00	Truck	Posture Constant Chair	900.00	UPS
Adjusting Cervical Chair	390.00	UPS	Manual Height Adjusting Stool	160.00	UPS
Knee Chest Table	700.00	UPS	Masters Spinalyzzer (Scolioscanner)	545.00	UPS 27.00
Gonstead Type Set	1550.00	Truck	Chiropractic Emblem-Aluminum	170.00	UPS 15.00
<i>(201 ABP, Chair & Knee Chest)</i>					

ACCESSORIES • CREPE HEADREST PAPER

DeJarnette Blockboard (16" X 20")	\$ 48.00	Child's Pelvic Blocks (Pair)	\$ 44.00
Large DeJarnette Blockboard (17" X 21")	53.00	Anterior Dorsal or Sacral Block	26.00
Pouch Blockboard	72.00	Anterior Dorsal Plate	38.00
Cranial Board	46.00	Low Silhouette Anterior Dorsal Plate	38.00
Sternal Roll (8 1/2" or 10 1/2")	23.00	Portable Face Plate (Open Bottom)	45.00
Footroll (Dutchman Roll)	32.00	Portable Face Plate (Solid Bottom)	38.00
Footroll (Flat Bottom)	35.00	Rectangular Pillow (15" or 18")	38.00
Adult Pelvic Blocks (Pair-Lloyd)	64.00	Vita Pedic Wedge	72.00
Adult Pelvic Blocks (Pair-DeJarnette Type)	62.00		

Minimum Shipping and Handling - \$10.00 - Per Order • Additional Charge Outside Contiguous U.S.A.

PORTABLE DROP HEADPIECES



SPECIFICATIONS:

	<u>Adult</u>	<u>Pediatric</u>
Base Width	5.5 or 9 Inches	6 Inches
Cushion Width	7 Inches	6 Inches
Length	9.5 Inches	6 Inches
Height	3.75 Inches	2.5 Inches

STANDARD FEATURES:

- Toggle Drop
- Solid Wood Construction
- Superior Crisp Drop
- Infinite Tension Control

PRICES: Pediatric	\$ 175.00
Adult - Narrow Base 5.5"	200.00
Adult - Wide Base 9	225.00
Extremity Drop - Deeper Drop w/ Contoured Cushion Top & Carry Handle.....	225.00

Narrow Base Adult Drop Head Shown

HEADREST PAPER - Hydrowoven (Tissue-Like • Soft • Strong)

One Carton (430 Feet Per Roll - 16 Rolls Per Carton)	\$ 85.00
Activator Headrest Paper (1,000 Sheets Per Box).....	50.00

Shipping & Handling Cost Included in Price of Paper. Additional Charge Outside Contiguous U. S.A.

EXAMINATION & TREATMENT TABLES



716	Examination & Therapy Table	\$ 545.00
716 P	Examination & Therapy Table	605.00
<i>716 & 716P - 28" Wide, 74" Long, 31" High - Shipping Weight 160 Lbs.</i>		
P - With Paper Roll & Dispenser		
P - Indicate direction of paper..Crosswise (Side to Side) or Lengthwise (Head to Foot)		

Items & Prices Subject to Change Without Notice Or Obligation

GENERAL INFORMATION

*Lloyd Table Company is proud of its record of high quality and dependable service.
Every order received, whether large or small, is given the same careful attention.*

NAUGA-PLUS, OUR STANDARD vinyl upholstery material is highly resistant to cracking from body oil and has a more pleasing touch to the skin. Cloth Upholstery also available for an additional charge.

STANDARD COLORS

Mahogany - CY-14
Chaps - CY 13
Buff - CY 15
Black - STA 22
Navy - STA 20
Spruce - CY-16

Gray - STA 12
Wine - STA 15
Dark Teal - STA 16
Black Cherry - STA 17
Raspberry - STA 18
Madder Brown - STA 21

SPECIAL ORDER COLORS

White - STA 11
Emerald - STA 13
Peacock - STA 14
Coastal Blue - STA 19
Cardinal - STA 23
Water - CY-09
Cabernet - PR 36

Laurel - PR 39

CLOTH COLORS

Mauve,Black,Navy,Gray
Teal,Tan,Burgundy,Dark
Brown
Dark Willow

NO REFUNDS OR CHANGES ON SPECIAL ORDER COLORS & CLOTH

FOAM: A variety of foam densities are available for doctor/patient preference. 3" of Foam standard (Galaxy, 400, Astro)
STANDARD FOAM DENSITIES: (As Used In The Past) Soft (35 Lb.); Medium (45 Lb.); Medium Firm (50 Lb.); Firm (70 Lb.); Extra Firm (90 Lb. Urethane)

SPECIAL ORDER FOAM DENSITIES: XXX Firm (Bonded); Maximum Firm (Closed Cell Foam) - *Call for Pricing - No Exchanges or Returns on Special Order Foam.*

* BIPA (Built in Paper Attachment)

International Surcharge on all electrically powered tables: \$200.00

- **PRIORITY RUSH: 10% ADDITIONAL CHARGE BASED ON TOTAL RETAIL COST OF ORDER**
The following products will ship within 12 working days - Galaxy & Astro Tables, Activator Type Elevation & Hylo, 400's, 401's & 402's, and Nicholas Side Posture. Does Not Include Special Order Items
All other products will ship within 4 working days. Does Not Include Special Order Items
- **NO PRIORITY RUSH ON GALAXY ULTIMATE TABLES.**
- **ALL PRIORITY RUSH ORDERS MUST BE PAID IN FULL PRIOR TO PROCESSING. ONCE AN ORDER HAS BEEN STARTED, ANY CHANGES, CANCELLATIONS, OR RETURNS WILL RESULT IN A FEE EQUAL TO THE PRIORITY RUSH CHARGE.**

- **NO CHANGES, RETURNS, OR REFUNDS ON SPECIAL ORDERED ITEMS.**
- **ANY CHANGES MADE TO A TABLE ORDER THAT HAS REACHED THE "WORK IN PROCESS" STAGE WILL BE CHARGED A MINIMUM OF \$100.00, NO MATTER WHAT THE CIRCUMSTANCES**
- **RESTOCKING CHARGE OF 10% ON ANY RETURNS (MINIMUM CHARGE OF \$10.00)**

CALL FOR A COPY OF OUR "RETURN POLICY"

If goods are damaged in transit, refuse to accept from the carrier until your delivery receipt is marked "Received in Damaged Condition". Lloyd Table Company will not be responsible for goods damaged in transit, all claims must be promptly filed with the transit company. If concealed damage is discovered upon unpacking, immediately notify the transportation company.

GUARANTEE

LLOYD equipment is guaranteed against defects in parts, materials (except upholstery) and workmanship as follows:

Galaxy Tables, Astro Tables, Activator Tables, 402 Tables, Cosmic III:

ONE YEAR (Full Parts)

SECOND & THIRD YEAR (Limited Parts)

Adjusting Bench Frames

LIFETIME

All Other Products:

THREE YEARS

Guarantee void if warranty card not returned to Lloyd Table Company